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Floatable special silicone for adapting of long-term noise and water plugs, addition curing


Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5


Fig. 6


Fig. 7

Fig. 8


Fig. 9


Fig. 10


Fig. 11

1. Mixing and Dosing

Extrude material from the cartridge with Automix1-system gun (Fig. 1). Place cartridge into mixing gun and remove closure cap by rotation. Eject a small amount of material (bleeding) until equal flow is being extruded from both orifices (Fig. 2). Mount the mixing cannula by using the guides on the cannula and the cartridge and lock in place by rotating in the opposite direction (Fig. 3). Feed material applying even pressure. Prior to use, dispense a small amount of material and visually check that the base and catalyst paste are mixed omogeneously (Fig. 4). Only then dosing can be performed individually. Leave mixing cannula on cartridge after use, thus sealing the cartridge.

2. Application using the direct method

Before taking an impression, the auditory canal and the eardrum have to be examined (Fig. 5). If distinctive features are found (e.g. inflammations or perforation of the eardrum, tubes) impression must not be taken. The auditory canal has to be cleaned and depilated. In case of insufficient cleaning, remaining cerumen can cause retarded setting resp. inhibition of setting process.

Before taking the impression, the eardrum must always be protected with an impression plug (e.g. blue secure) placed at the end of the external auditory canal in front of the eardrum (Fig. 6). Now **aquaplust** is injected directly into the prepared ear of the patient (Fig. 7) as described under point 1. Keep the end of the cannula immersed on the surface of the material during dispensing. After curing remove the impression carefully from the ear (Fig. 8). Finally it is essential to re-examine the ear to make sure it has not been damaged (Fig. 9). The earmould is trimmed and polished as described under point 4.

3. Indirect application in the laboratory

Prepare the ear impression for the post-processing, apply wax or lacquer. All commercially available plasters and gel materials can be used for the fabrication of the negative form. Coat the plaster negative form with a usual alginate based insulation. An insulation of gel forms is not necessary. The bubblefree injected material vulcanises at room temperature without pressure application. If a handle is needed, we recommend the installation of **detax handy** or **grip tool** according to special instructions. A plaster negative form with counter can be made for the fabrication of otoplastics, in order to shorten trimming: Fill one part of the flask with plaster and embed the insulated impression. After setting of the plaster open the flask. Scald both parts of the flask with boiling water and coat them with a usual alginate based insulation. Extrude **aquaplust** from the cartridge (described under point 1) into the negative form and close flask with counter accurately (Fig. 10). After vulcanisation open the flask and remove the die.

4. Final trimming and varnishing

For moulding and surface trimming use special cutters or grinding sleeves. If a handle is needed, we recommend the installation of **detax handy** or **grip tool** according to special instructions. For smoothing the surface, use the **DETAX silicone lacquers** according to their special instructions (Fig. 11).

Important working hints

- Do not combine with condensation curing silicones.
- Cured silicone materials are chemically inert – spots on clothing should be avoided.
- Where movement is impeded or blocked, discard the cartridge. Do not apply force.
- Be sure to remove all traces of silicone material from the patient's ear.
- Latex gloves and latex contaminated surfaces as well as cerumen, cremes and resins may inhibit the setting reaction of **aquaplust** (we recommend standard gloves made of nitrile rubber or polyethylene).
- A careful preparation (point 1) as well as a thorough cleaning of the auditory canal is required before taking an impression.
- After each impression taking the auditory canal has to be re-examined!
- Please follow the instructions of the safety data sheet!

Safety advice

The impression taking of the external auditory canal has to be executed by trained specialist staff only. The working instructions and precautionary measures have to be strictly observed. Non-observance might lead to irreparable damage of the ear or ear drum. DETAX is not liable for any damage caused by improper application of the impression material.

Further information:

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

For use by trained specialists only.

Indications for use:

- Direct or indirect fabrication of noise and water protection earmoulds

Technical data:

- **Mixing volume:** 50 ml (cartridge)
- **Mixing ratio:** 1:1
- **Colour code:** base: red, blue, pink, orange, yellow, green catalyst: white
- **Mixing time:** omitted (Automix1-system)
- **Setting time:** **direct method:** approx. 5 min. at body temperature **indirect method:** approx. 10 - 15 min* at room temperature (approx. 23 °C / 73 °F)
- **Final hardness:** approx. 35 Shore A
- **Strain in compression:** approx. 7 %
- **Recovery from deformation:** approx. 99.5 %
- **Linear dimensional change:** approx. 0.3 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**

Ordering information:

aquaplust
Standard packing 02423

8 x 50 ml cartridges
base + catalyst, colour
assorted: 4 x red, 4 x blue

aquaplust colormix
Standard packing 03205

8 x 50 ml cartridges
base + catalyst
assorted: 2 x red, 2 x blue,
2 x orange, 1 x yellow,
1 x green
32 mixing cannulas optima

Mixing cannulas optima

6 mm, 50 pcs. 03187
6 mm, 100 pcs. 03188

express coat 03012
Aluminium bottle, 100 ml

dip coat 03014
Aluminium bottle, 100 ml

super coat 02949
Aluminium bottle, 100 ml

aqua coat 03740
Aluminium bottle, 2 x 50 ml

micropor lacquer 02169
aluminium bottle, 250 ml

detax handy 02469
5 x blue, 5 x red

grip tool 03083
2 x blue, 2 x red

impression plugs 03321
blue secure
100 pcs., conical

DETAX
GmbH & Co. KG
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
Telefon: 0 72 43/510-0 · Fax: 0 72 43/510-100
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a licensed healthcare practitioner (or trained specialist personnel).

Made in Germany

01/2019